


WORTH TRUST - ASSISTIVE DEVICES FOR THE VISUALLY IMPAIRED


| Code No. | Product | Description | Image |
|----------|-------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| B1 | Perkins Brailler - Standard Classic Model | Perkins Brailler can emboss 25 lines with 42 cells on 11" x 8 ½ " sheet of paper. Includes leather dust cover and wooden eraser. Manually operated braille teaching and writing device. Highly durable. Students can take notes or exams quickly in order keep up with sighted peers, increasing their motivation and confidence. |  |
| B2 | Perkins Brailler - SMART Model | The Perkins SMART Brailler makes learning braille fun and easy. Teachers and parents can observe letters on the display screen, while the student is motivated by the letter or word being spoken and produced in braille. A video screen that displays SimBraille and large print, combined with audio feedback, for an instantaneous, multi-sensory learning experience. |  |
| B3 | Perkins Brailler - Large Cell Model | Easier to Read: Produces large braille cells with more space between them so that people with tactile problems can more easily read and produce braille. |  |
| B4 | Perkins Brailler - Unimanual Model | Designed for One Hand: The Perkins Unimanual Brailler works like other braille writers, except the keys for dots 1,2 and 3 will remain in a locked position until dots 4,5 and 6 are added or until the spacebar is depressed. |  |



| | | | |
|------|-----------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| B5 | Perkins Brailler - Electric Model | Easy Writing: A lighter touch goes a long way with the Electric Blue Perkins Brailler. Students can write longer and faster. It's also great for children as well as adults with weaker hands. Felt pad, dust cover and wooden eraser included. UL and CSA approved |  |
| BC1 | Brailler Leather Dust Cover | This black leather embossed dust cover helps protect your Perkins Brailler from dirt and dust. The Brailler handle is accessible through the slot in the dust cover, promoting portability. Cover fits all Perkins Brailler. |  |
| BC2 | Brailler Soft Side Carry Case | The Soft side Carrying Case has two key-operated locks to keep the brailler secure during transit. It is made from a lightweight, water-resistant material with a dual, heavyweight zipper. |  |
| UBK1 | Universal Braille Kit - Standard Model | Braille Slate A4 size with one Stylus, English & Hindi / Tamil / Telugu Alphabet Trainer Plate, Geometry Set Large, White Cane - 4 folding, Braille Word Forming Block (5 letters), Braille Cube, Tracer Wheel, Braille Paper (20 Sheets), Taylor's Mathematical Frame in plastic with Type sets in plastic, Braille Abacus, Manual on how to use each of the assistive devices, Back pack. |  |
| UBK2 | Mini Universal Braille Kit for Beginner | Braille Slate A4 size with one Stylus, Geometry Set Large, Vikas Composer, Braille Word Forming Block (5 letters), Braille Cube, Tracer Wheel, Braille Paper (20 Sheets), Mini Taylor's Mathematical Frame in plastic with cubes in plastic, Braille Abacus, Manual on how to use each of the assistive devices, Bag. |  |



| | | | |
|----|------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| S1 | Stylus - Flat Type - Plastic | Write Braille in any Language - 70 mm in length |  |
| S2 | Stylus - Knob Type - Plastic | Write Braille in any Language - 60 mm in length |  |
| S3 | Stylus - Saddle Type - Plastic | Write Braille in any Language - 60 mm in length |  |
| S4 | Stylus - Reversible Type - Plastic | Write Braille in any Language and also has an eraser |  |
| S5 | Braille Slate A4 size - Grouse model - Plastic | 30 cells, 27 lines. This is a Braille Slate meant for Beginners. A thick paper sheet of the same size as the slate should be inserted between the two leaves of the slate. Braille text in any language should then be embossed on the paper sheet, one letter in each cell of the slate using the stylus. Each cell has three notches on each of its vertical edges, to control the exact position of the stylus while embossing the dots in corresponding positions. It must be remembered that the embossing of text has to be from right to left starting from the right edge of the sheet. |  |


| | | | |
|-----|-----------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| S6 | Braille Slate A4 size - Standard Plain model | 30 cells, 27 lines. A thick paper sheet of the same size as the slate should be inserted between the two leaves of the slate. Braille text in any language should then be embossed on the paper sheet, one letter in each cell of the slate using the stylus. It must be remembered that the embossing of text has to be from right to left starting from the right edge of the sheet, with the individual letter codes being embossed laterally inverted. This is necessary as the paper has to be turned over for the embossed matter to be read by touch. |  |
| S7 | Braille Slate Pocket size - Plastic | 21 cells, 9 lines. 10 cms. x 17 cms. A thick paper sheet of the same size as the slate should be inserted between the two leaves of the slate. Braille text in any language should then be embossed on the paper sheet, one letter in each cell of the slate using the stylus. It must be remembered that the embossing of text has to be from right to left starting from the right edge of the sheet, with the individual letter codes being embossed laterally inverted. This is necessary as the paper has to be turned over for the embossed matter to be read by touch. |  |
| E2 | Eraser made out of Plastic | Quick Fix: Fix braille mistakes fast with the Eraser. Use the blunt tip to rub out unwanted dots. No need to start over. |  |
| BP1 | Braille Paper 100 sheets, A4 size, White (140 GSM) | Used for writing in a Braille Slate and Perkins Brailler |  |
| BP2 | Braille Paper 100 sheets, 11 x 11 ½ size, White (140 GSM) | Used for writing in a Braille Slate and Perkins Brailler |  |

| | | | |
|-----|--------------------------------------------------|-------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| TA1 | Braille Letter Cube - Plastic | Rotate the cube's segments to form different Braille letters |  |
| TA2 | Braille Word Forming Block - 5 letters - Plastic | To form Braille words upto 5 letters in any language. |  |
| TA2 | Braille Word Forming Block - 8 letters - Plastic | To form Braille words upto 8 letters in any language. |  |
| G1 | Geometry set - Standard in box - Plastic | 30° & 45° Set squares, Compass, Scale-15 Cms. & Protractor in a box |  |
| G2 | Geometry set - Large - Plastic | 30° & 45° Set squares, Compass, Scale-30 Cms. & Protractor in a plastic cover |  |

| | | | |
|-----|-------------------------------------------|---------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| G4 | Braille Scale 12" | Made of plastic - 12" |  |
| SW1 | Tracing Wheel - Plastic | It is a drawing tool to make embossed dotted lines on Plain paper which can be sensed by touch. |  |
| RM1 | Rubber Mat - 5 mm | Writing and Drawing surface. |  |
| A1 | Abacus - Pocket size - Plastic - with box | 15 vertical bars - Frame 90 mm x 160 mm packed in Plastic Box |  |
| M1 | Taylor's Mathematical Frame - Plastic | 18 lines x 22 cells in each line. It provides a complete tactile solution for the Visually Impaired to learn mathematics. |  |

| | | | |
|------|-------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| M2 | Type sets for Mathematical Frame, Plastic | 100 pieces packed in a plastic box. Both Arithmetic and algebraic types to be used along with the Taylor Frame. |  |
| M3 | Mini Taylor's Mathematical Frame - Plastic with cubes | Used by the beginners to learn mathematics. |  |
| TA 4 | English Alphabet Trainer Plate - Plastic | To learn English alphabets and numerals by writing over the space provided on the plate |  |
| TA5 | Tamil Alphabet Trainer Plate - Plastic | To learn Tamil alphabets and numerals by writing over the space provided on the plate |  |
| TA6 | Hindi Alphabet Trainer Plate - Plastic | To learn Hindi alphabets and numerals by writing over the space provided on the plate |  |

| | | | |
|-----|-----------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| TA7 | Telugu Alphabet Trainer Plate - Plastic | To learn Telugu alphabets and numerals by writing over the space provided on the plate |  |
| TA8 | Arabic Alphabet Trainer Plate - Plastic | To learn Arabic alphabets and numerals by writing over the space provided on the plate |  |
| WB1 | Word Building Educational Kit - Plastic | A board of plastic with 60 pieces of alphabets in Braille and in English that can be interlocked to form words. |  |
| V1 | Vikas Composer | To form braille sentences in any language. |  |
| D1 | Dominoes Double Six | Tactile Dominoes Double Blank. A standard 28-piece set with the highest being the Double Six, made of plastic. The black dots are raised for easy identification of each piece. Each domino measures 2.5 cm x 4.7 cm. |  |

| | | | |
|------|------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| SG1 | Signature Guide - Aluminium | To help the Visually Impaired to sign in the specified space. |  |
| CRS1 | White Mobility Cane - Rigid | Made of aluminium powder coated with plastic tip and coloured reflective marking in red. Length: 110, 120, 130 Cms or as per customer requirement. Handle - with or without rubber handle. |  |
| CFS1 | White Mobility Cane - 4 folding - Powder coated | Made of aluminium powder coated with plastic tip and coloured reflective marking in red. Length: 110, 120, 130 Cms or as per customer requirement. Can be folded and carried in a bag. |  |
| CFS2 | White Mobility Cane - 4 folding - Anodized | Made of aluminium anodized with plastic tip and coloured reflective marking in red. Length: 110, 120, 130 Cms or as per customer requirement. Can be folded and carried in a bag. |  |
| CFS3 | White Mobility Cane - 4 folding - Premium Model - Anodized | Made of premium aluminium anodized with plastic tip and coloured reflective marking in red. Length: 110, 120, 130 Cms or as per customer requirement. Can be folded and carried in a bag. |  |

| | | | |
|------|-------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| CFS4 | White Mobility Cane - 5 folding - Powder coated | Made of aluminium power coated with plastic tip and coloured reflective marking in red. Length: 110, 120, 130 Cms or as per customer requirement. Can be folded and carried in a bag. |  |
| RS1 | Low Vision Reading Stand | Made of Metal and Plastic Foldable. Size 35 Cms x 25 Cms. It can be adjusted to the required angle for easier reading of persons with low vision. | |

For Ordering please contact:

Nunes Ignatius Myke
Manager - Brailler Division & Marketing Services
WORTH Trust
48, New Thiruvalam Road, Katpadi, Vellore - 632007, Tamil Nadu,
India.
Mobile: 9442592306
Tel: +91 416 2242739
Mail: myke@worthtrust.org.in
Web: www.worthtrust.in